


MyDay[®] daily disposable toric



High-performing 1 Day contact lenses
for wearers with astigmatism.



You demand more from your toric contact lenses: Beyond clear, stable vision, you insist on extraordinary comfort to help you stay on top of your high intensity lifestyle. MyDay[®] daily disposable toric lenses from CooperVision, a world leader in toric contact lenses, are made with the same design features found in Biofinity[®] toric. They are also part of the MyDay daily disposable family: Our softest silicone hydrogel lenses ever* and the world's first 1 Day contact lenses made with CooperVision[®] Aquaform[®] Technology to help keep your eyes looking clear and white and feeling comfortable all day long. With MyDay daily disposable toric, you're ready to take on your day.

Replacement Schedule



Corrects



Features & Benefits

- Excellent vision correction and a proven, stable and comfortable fit if you have astigmatism. CooperVision[®] toric contact lenses benefit from our expertise in toric design with Optimized Toric Lens Geometry™
- Long-lasting comfort: advanced Aquaform[®] Technology locks in water—so you can have incredible long-lasting comfort.
- High breathability: silicone hydrogel material allows 100% oxygen to reach your eyes¹, to help keep them clear and white.**
- UV blocking: helps protect your eyes from ultraviolet rays, helping to support eye health.†

*Compared among CooperVision[®] 1 day silicone hydrogel contact lenses. Data on file.

**High-oxygen transmissibility promotes clear, white eyes.

†Warning: UV blocking contact lenses are not substitutes for protective UV-absorbing eyewear, such as UV-absorbing goggles or sunglasses, because they do not completely cover the eye and surrounding area. Patients should continue to use UV-absorbing eyewear as directed

References:

1. Brennan, NA. Beyond Flux: Total Corneal Oxygen Consumption as an Index of Corneal Oxygenation During Contact Lens Wear. Optom Vis Sci. 2005;82(6):467-472.

©2021 CooperVision

MyDay[®], CooperVision[®], and Aquaform[®] Technology are registered trademarks of The Cooper Companies, Inc., and its subsidiaries.